
SAXS Documentation

Release 2.0

Christian Meisenbichler


May 12, 2015

CONTENTS

1	Install	3
2	The Tools	5
2.1	The Saxsdog	5
2.2	The Converter	5
2.3	Plotchi	6
2.4	Saxsdmerge	6
3	The Saxsdog Network	9
3.1	The SAXSNetwork configuration	9
3.2	The Saxsdog Server	10
3.3	The Saxs Leash	11
3.4	Saxs Leash Commandline	11
3.5	The Saxsdog Network Protocol	13
3.6	The Protocol Schemas	17
4	The Dedector Calibration File	33
4.1	Calibration file Reference	33
5	The Technology	43
5.1	Integration as Matrix-Vector Multiplication	43
5.2	Over sampling	44
5.3	The Geometry	44
5.4	Polarization Correction	48
5.5	Compare With Fit2d	51
5.6	Integrating a Constant Image With Masked Values	55
5.7	Statistics	58
5.8	GISAXS Slices	60
6	SAXSDOG Developer Documentation	61
6.1	The Leash	61
6.2	The Image Queue	64
6.3	Data Consolidation or Datamerger	64
7	SAXS Module API	65
	Index	69



The SAXS Python package implements analysis tools for Small Angle X-Ray Scattering (SAXS) data. The first and most important one is to efficiently integrate 2d sensor data to an angle dependent diffraction curve.



The SAXS module consists of a Python library and 3 command line tools: *The Saxsdog*, *Plotchi* and *The Converter* and the *The Saxsdog Network* that integrates all of it and provides a GUI, *The Saxs Leash*.

INSTALL

The SAXS Package is distributed as a Python package. So in order to use it, you need a Python system installed. It depends on following Python modules that don't come with the standard Python:

```
numpy scipy matplotlib jsonschema bitarray watchdog sphinxcontrib-programoutput\  
sphinxcontrib-programsscreenshot pyzmq prettyplotlib pandas
```

They are all available through “pip”, so, the command:

```
>>pip numpy scipy matplotlib jsonschema bitarray watchdog sphinxcontrib-programoutput\  
sphinxcontrib-programsscreenshot pyzmq prettyplotlib pandas
```

should get all the modules. For Windows, use the Anaconda Python distribution which includes pip.

The code can be obtained on github: <https://github.com/ChristianMeisenbichler/SAXS> where you would also find a “Download Zip” button. After unpacking or cloning with git you end up with a directory called “SAXS” containing the files. Go there, and type into the command line:

```
python setup.py install
```

This installs the Python module to the environment, creates the command line tools and installs them on the system. Where that is, depends on the Python installation. The setup script will also try to satisfy all the dependencies by downloading and installing the missing packages.

2.1 The Saxsdog

The saxsdog is a script that converts directories with images to curves. It can use multiple threads and watch the file system for changes.

For help on the usage type:

```
$ saxsdog --help
Usage: saxsdog [options] directory/to/watch
```

Options:

```
-h, --help                Show this help message and exit.
-c FILE, --calibration=FILE Path to calibration file (JSON).
-t THREADS, --threads=THREADS
                          Number of concurrent threads.
-m, --plotmonitor         Show a live updating plot window.
-w, --watch                Watch directory for changes, using file system events
                          recursively for all sub directories.
-r, --resume              Skip files that are already converted.
-o OUTDIR, --out=OUTDIR   Specify output directory.
-R RELPATH, --relpath=RELPATH
                          Specify output directory as relative path to image
                          file. Default: '../work'
-s, --svg                  Write plot to svg file.
-p, --png                  Write png of original.
-P, --profile              Make a time Profile and print it.
-S, --silent               Less output.
-n, --nowalk               Don't scan for files already there, only watch file
                          system if -w flag is given.
-D, --walkdirinthreads    Search all directories in parallel process.
-V, --servermode           Servermode.
```

The calibration file must be a valid *Calibration file Reference*

2.2 The Converter

The converter extracts information from the calibration.txt generated by fit2d and adds them to a SAXS.calibration configuration file. (*Calibration file Reference*)

```
$ saxsconverter --help
Usage:
To extract data from Fit2d output:
  saxsconverter [options] calibration.txt ouput.saxsconf
```

Or to convert fom older saxsconf:

```
saxsconverter [options] cal.saxsconf ouput.saxsconf
```

Options:

```
-h, --help          show this help message and exit
-t FILE, --template=FILE
                    Path to calibration file which serves as template.
```

If there is a target file and it is a valid *Calibration file Reference*, the parsed values are added or updated in place.

2.3 Plotchi

The tool “plotchi” plots a list of “.chi”-files:

```
$ plotchi --help
```

```
Usage: plotchi [options] CHIFILE [List of more ".chi" files]
```

Options:

```
-h, --help          show this help message and exit
-o FILE, --out=FILE Write the plot to FILE. The format is derived from the
                    suffix, e.g. '.svg', '.pdf'.
-c, --compare       Compare datasets to first one.
-l, --log           Use log scale.
-n, --no-legend     Hide legend.
-t TITLE, --title=TITLE
                    Give plot title.
-s N, --skip=N      Skip first N points.
-k N, --clip=N      Clip last N points.
-x TYPE, --xaxisistype=TYPE
                    Select type of X axis scale, might be
                    [linear|log|symlog]
-y TYPE, --yaxisistype=TYPE
                    Select type of Y axis scale, might be
                    [linear|log|symlog]
```

2.4 Saxsdmerge

This tool addresses a very specific problem only applicable under our beamlines needs.

```
$ saxsdmerge --help
```

```
Usage: saxsdmerge [options] iMPicture/dir peakinteg.log datalogger.log
```

Options:

```
-h, --help          show this help message and exit
-t SEC, --timeoffset=SEC
                    Time offset between logging time and time in
                    imagedata.
-l, --syncfirst     Sync time by taking the time difference between first
                    shutter action and first image.
-o FILE, --outfile=FILE
                    Write merged dataset to this file. Format is derived
                    from the extesion.(.csv|.json|.hdf)
-b, --batch         Batch mode (no plot).
-c, --includechi    Include radial intensity data (.chi) in hdf.
-f, --includetif    Include all image data in hdf.
-C FILE, --conf=FILE
                    Use config in FILE to merge the data (ignore other
                    options)
```


What it does is the following: It merges the two log files by interpolating the missing data in the joined table of the the 1 sec interval data logger and the shutter log. The interpolation method takes the closest previous entry. The two logfiles are presumed to have the same clock. This data is then merged with the data extracted from the image headers and dedector logs and is written to a output file as a Table. As the image time stamps are from a different clock which might have a significant offset, the `-l` and `-t` options allow for dealing with that. In order to check if the time synchronization is reasonable, the tool displays a graph with the shutter times and the exposure times from the images.

One of the output file format options is hdf, the hirarchical data format, which ia a way to create portable binary data files. If you specify hdf, the radial average function `".chi"` and even the image data can be written to the hdf file. The hdf export is done with pandas hdf exporter which in turn uses `PyTables`. The `".chi"` files reside under the "Graphs" Label, the log data under the "LogData" label and the images if present under the "Images" Label. It can be restored to a Pandas pannel or data frame with:

```
import pandas as pd
store=pd.HDFStore("/home/chm/saxs/merged.hdf")
Curves=store["Curves"] # gives Pandas panel with all the curves
Data=store["Data"] # gives Pandas data frame of merged log and image parameters
```


THE SAXSDOG NETWORK

The network may consist of 3 different services. The “Saxsdog Server” does the image processing. The “Saxs Feeder” publishes new file events and the “Saxs Leash” controls and configures the server. The feeder is optional.



3.1 The SAXSNetwork configuration

The Saxsdog server and the Saxsleash have a common configuration file, which tells them how to connect with each other and which also includes a shared secret for authentication. If you want two computers to connect via the Saxsleash you need to have a copy of the file on each of them.

To create such a configuration, use the command:

```
$ saxsnetconf
```

It will ask for the Feeder URL and for the Saxsdog Server URL. Then it will generate a random secret and save the file in file in `$Home/.saxsdognetwork`. You will have to copy the file to the other computers you need to allow to connect to your network. The secret must be the same on all of them.

```
[{
  "Server": "tcp://hostname:port",
  "Feeder": "tcp://hostname:port",
  "Secret": "Some large random string."
  "Name": "Pilatus1M"
}]
```

The authentication is done by hashing the request and the secret including a time stamp. The time stamp is checked if it lies within 900 seconds of the servers time.

3.2 The Saxsdog Server


The Saxdog Server is the program that is started on the processing computer (node). It may subscribe to a “new file” event service.

```
$ saxsdogserver --help
Usage: saxsdogserver [options] basedir
```

Options:

```
-h, --help          show this help message and exit
-p port, --port=port  Port to offer command service. Default is 7777.
-t THREADS, --threads=THREADS
                    Number of concurrent processes.
-f tcp://hostname:port, --feeder=tcp://hostname:port
                    Specify the URL of the new file event service (Saxsdog
                    Feeder)
-w, --watch          Watch directory for changes, using file system events
                    recursively for all sub directories.
-R RELPATH, --reldir=RELPATH
                    Specify output directory as relative path to image
                    file. Default: '../work'
-o OUTDIR, --out=OUTDIR
                    Specify output directory
-d, --daemon         Start server as daemon
```

3.3 The Saxes Leash



The Leash Program is a GUI to load calibrations into the Saxdog Server and monitor the processing of the data. It provides a calibration editor, as mask preview and basic data import from *saxsconverter*.

The main window has 3 tab cards. The first is for setting up the server, the second to review the currently processed data and the third for basic statistics. The command to launch it is.

```
leash
```

3.4 Saxes Leash Commandline

The “Saxes Leash” command line client can issue the commands for the Saxsdog Server.

```
$ saxsleash --help
```

```
Usage: saxsleash close|abort|new|get|plot|plotdata|readdir|stat|listdir|putplotdata|fileslist|mer
```

Options:

```
-h, --help          show this help message and exit
-S tcp://HOSTNAME:PORT, --server=tcp://HOSTNAME:PORT
                    URL of "Saxsdog Server"
-s N, --skip=N      plot: Skip first N points.
-k N, --clip=N     plot: Clip last N points.
-x TYPE, --xaxisistype=TYPE
                    plot: Select type of X axis scale, might be
```

```
                                [linear|log|symlog]
-y TYPE, --yaxsisstype=TYPE      plot: Select type of Y axis scale, might be
                                [linear|log|symlog]
-N SERVERNO, --serverno=SERVERNO
                                select server from config list by index default:0
```

Most of the command line options are about the `plot` command, but in order to visualize the processed data, one has to send the commands to setup a calibration.

3.4.1 New

```
$ saxsleash new cal.json data/AAA_integ.msk data/
```

The `new` command loads a calibration and starts the queue to receive new files. It requires 3 arguments:

1. Calibration file. as in *The Dedector Calibration File*,
2. mask file,
3. directory where the image files are or are going to be.

If there is a queue running, this command will abort the other one and replace it. One server can have only one queue at a time.

3.4.2 Plot

```
$ saxsleash plot
```

The `plot` command will grab the next image and show a plot of the result in a window. This command will be repeated until the user interrupts it with `Ctrl-C`.

3.4.3 Close

```
$ saxsleash close
```

Closes the queue. Which means, the server will process what is left in the queue but ignore all new files.

3.4.4 Abort

```
$ saxsleash abort
```

The `abort` command will close the queue and stop all data processing processes. It will only wait for each process to finish the picture they started before. The remaining pictures in the queue are ignored.

3.4.5 Stat

```
$ saxsleash stat
```

Return basic statistics data about the processes.

3.4.6 Read Dir

```
$ saxsleash readdir
```

This command will put all the images in the configured directory into the queue. This is useful to reprocess pictures.

3.4.7 List Directory

```
$ saxsleash listdir Serverdirectory
```

Returns a list of files and directories on the server. Mostly used for the Leash GUI dialogues.

3.4.8 Files List

```
$ saxsleash fileslist
```

Returns a list of processed files.

3.4.9 Merge Data

```
$ saxsleash mergedata mergeconf.json
```

Sends command to merge data. Takes a merge configuration file as argument.

3.5 The Saxsdog Network Protocol

3.5.1 The Saxsdog Server

The saxdog server can watch for files system events for himself or subscribe to a zmq service, the Saxsdog Feeder, that publishes new file names. The server can process the new images according to one calibration. The server may only have one calibration at a time, it is not designed to be used by multiple users at the same time.

3.5.2 The Saxsdog Feeder

The “Saxsdog Feeder” service offers file events for subscription. It should not do any buffering or pre-selection, just send a new message when any new file was copied and is ready for processing. Also when a file is overwritten: Send a message. It should however, only send this event, when the file is completely written to the file system.

New file events are composed of the following message:

```
{
  "command": "new file",
  "argument": "/Path/to/file/"
}
```

The service must be a ZeroMQ zmq.PUP socket. This code is a simulation of the messages:

```
import zmq
import random
import sys
import time
import os
import json
from optparse import OptionParser
```

```
def startfeeder():
    """
    Simulator for new file announcement service. For development and testing.
    """
    parser = OptionParser()
    usage = "usage: %prog [options] "
    parser = OptionParser(usage)
    parser.add_option("-p", "--port", dest="port",
                    help="Port to offer file changes service", metavar="port", default="")
    parser.add_option("-d", "--dir", dest="dir",
                    help="Directory to monitor", metavar="dir", default=".")
    parser.add_option("-s", "--sdir", dest="sdir",
                    help="server dir, (prefix to filepaths)", metavar="dir", default=".")
    (options, args) = parser.parse_args(args=None, values=None)

    context = zmq.Context()
    socket = context.socket(zmq.PUB)
    if options.port=="":
        conf=json.load(open(os.path.expanduser("~/+os.sep+".saxsdognetwork")))
        port=conf['Feeder'].split(':')[1]
    else:
        port=options.port
    print "connecting:", "tcp://*:%s" % port
    socket.bind("tcp://*:%s" % port)

    fileslist=[]
    if len(args)>0 and options.dir=="." :
        dirtosearch =args[0]
    else:
        dirtosearch =options.dir
    for path, subdirs, files in os.walk(dirtosearch):
        for name in files:
            if name.endswith('tif'):
                fileslist.append( os.path.join(path, name))
    messageobj={"command":"new file","argument":""}
    while True:
        for file in fileslist:
            print file
            messageobj['argument']=file
            message=json.dumps(messageobj)
            socket.send(message)
            time.sleep(7)

if __name__ == '__main__':
    startfeeder()
```

3.5.3 The Saxsdog Leash

The Saxsdog Leash is a user-facing control interface. There, the user should enter new calibrations and specify the data directories connected to it. During the processing, it shows a graph of one of the current images.

It may send the following commands:

Close

Request:

```
{
  "command":"close queue",
  "time":1404979588.715198,
  "sign":"Signature generated for request"
}
```

Answer:

```
{
  "result":"queue closed",
  "data":{
    "stat": {
      "time interval": 0.8776118755340576,
      "queue length": 0,
      "frames per sec": 10.25510279760422,
      "images processed": 235, "pics": 9
    }
  }
}
```

Abort

Request:

```
{
  "command":"abort queue",
  "time":1404979588.715198,
  "sign":"Signature generated for request"
}
```

Answer:

```
{
  "result":"queue stopped emptied and closed",
  "data":{
    "stat": {
      "time interval": 0.8776118755340576,
      "queue length": 0,
      "frames per sec": 10.25510279760422,
      "images processed": 235,
      "pics": 9
    }
  }
}
```

New

Request:

```
{
  "command":"new queue",
  "argument":{
    "directory":["path","to","data"],
    "calibration":{},
    "maskbin":""
  },
  "time":1404979588.715198,
  "sign":"Signature generated for request"
}
```

Answer:

```
{ "result":"new queue",
  "data":{
    }
}
```

Plot

Request:

```
{ "command":"send plot",
  "time":1404979588.715198,
  "sign":"Signature generated for request"
}
```

Answer:

```
{
  "result":"plot data",
  "data":{
    "filename":"/name/.tiv" ,
    "stat": {
      "time interval": 0.8776118755340576,
      "queue length": 0,
      "frames per sec": 10.25510279760422,
      "images processed": 235,
      "pics": 9
    },
    "array":[[0],[0],[0]]
  }
}
```

Readdir

This puts all existing files in the queue directory into the queue again.

Request:

```
{
  "command":"readdir",
  "time":1404979588.715198,
  "sign":"Signature generated for request"
}
```

Answer:

```
{
  "result":"directory refilled queue",
  "data":{
    "stat": {
      "time interval": 0.8776118755340576,
      "queue length": 0,
      "frames per sec": 10.25510279760422,
      "images processed": 235, "pics": 9
    }
  }
}
```

Stat

Get basic processing statistics.

Request:

```
{ "command":"stat", "argument":{},
  "time":1404979588.715198,
  "sign":"Signature generated for request"}
```

Answer:

```
{
  "data": {
    "stat": {
      "time interval": 711.6886098384857,
      "queue length": 0,
      "frames per sec": 9.972057866165134,
      "images processed": 7332,
      "pics": 7097
    }
  },
  "result": "stat"
}
```

Error

In case of error in the Saxsdog server it will return an error message:

```
{
  "result":"Error",
  "data":{"Error":"Error message"}
}
```

3.6 The Protocol Schemas

3.6.1 Leash Request Schema

This describes how the requests to the server must be composed.

The '*' signifies a required Field.

Schema for requests from Saxes Leash to Saxes Server

Type object

Contains *command**, *argument*, *sign*, *time*

Required True

JSON Path

• #

Example JSON:

```
{"command": "close"}
```

command

Type string

values [close, abort, new, get, plot, plotdata, readdir, stat, listdir, fileslist, mergedata, getmergedata]

Required True

JSON Path

- # ['command']

Example JSON:

```
{"command": "close"}
```

argument

Type object

Contains *calibration, mergeconf, data, directory*

Required False

JSON Path

- # ['argument']

Example JSON:

```
{"argument": {}}
```

calibration

Calibration data according to *The Dedector Calibration File*

Type object

Contains /

Required False

JSON Path

- # ['argument']['calibration']

Example JSON:

```
{"calibration": {}}
```

mergeconf

Datamerger Configuration

Type object

Contains /

Required False

JSON Path

- # ['argument']['mergeconf']

Example JSON:

```
{"mergeconf": {}}
```

data

type object

Required False

JSON Path

- # ['argument']['data']

Example JSON:

```
{"data": {}}
```

directory

type object

Required False

JSON Path

- # ['argument']['directory']

Example JSON:

```
{"directory": {}}
```

sign

Signature of request

Type string

Required False

JSON Path

- # ['sign']

Example JSON:

```
{"sign": ""}
```

time

time in seconds (pythons time.time())

Type number

Required False

JSON Path

- # ['time']

Example JSON:

```
{"time": 0}
```

3.6.2 Leash Response Schema

This is the Schema of the data the server sends back as a response.

The '*' signifies a required Field.

Schema for requests from Saxes Leash to Saxes Server

Type object

Contains *result**, *data**

Required True

JSON Path

- #

Example JSON:

```
{"data": {}, "result": ""}
```

result

Type string

Required True

JSON Path

- # ['result']

Example JSON:

```
{"result": ""}
```

data

Type object

Contains *cal*, *Error*, *synplot*, *directory*, *attachments*, *threads*, *dircontent*, *history*, *fileslist*, *stat*, *filename*, *IntegralParameters*, *graphs*

Required True

JSON Path

- # ['data']

Example JSON:

```
{"data": {}}
```

cal

Type object

Required False

JSON Path

- # ['data']['cal']

Example JSON:

```
{"cal": null}
```


Error

type object

Required False

JSON Path

- # ['data']['Error']

Example JSON:

```
{"Error": {}}
```

syncplot

type object

Required False

JSON Path

- # ['data']['syncplot']

Example JSON:

```
{"syncplot": {}}
```

directory

Directory this queue is going to use. New files in other directories are going to be ignored.

Type array() items:

Required False

Default [u'.', u'', u'']

JSON Path

- # ['data']['directory']

Example JSON:

```
{"directory": [".", "", ""]}
```

attachments

type object

Required False

JSON Path

- # ['data']['attachments']

Example JSON:

```
{"attachments": {}}
```

threads

Type integer

Required False

JSON Path

- # ['data']['threads']

Example JSON:

```
{"threads": 0}
```

dircontent

type object

Required False

JSON Path

- # ['data']['dircontent']

Example JSON:

```
{"dircontent": {}}
```

history

type object

Required False

JSON Path

- # ['data']['history']

Example JSON:

```
{"history": {}}
```

fileslist

type object

Required False

JSON Path

- # ['data']['fileslist']

Example JSON:

```
{"fileslist": {}}
```

stat

type object

Contains *queue length, images processed, time, start time, mergecount*

Required False

JSON Path

- # ['data']['stat']

Example JSON:

```
{"stat": {}}
```

queue length

Type integer

Required False

JSON Path

- # ['data']['stat']['queue length']

Example JSON:

```
{"queue length": 0}
```

images processed

Type integer

Required False

JSON Path

- # ['data']['stat']['images processed']

Example JSON:

```
{"images processed": 0}
```

time

Type number

Required False

JSON Path

- # ['data']['stat']['time']

Example JSON:

```
{"time": 0}
```

start time

Type number

Required False

JSON Path

- # ['data']['stat']['start time']

Example JSON:

```
{"start time": 0}
```

mergecount

Type number

Required False

JSON Path

- # ['data']['stat']['mergecount']

Example JSON:

```
{"mergecount": 0}
```

filename

Type string

Required False

JSON Path

- # ['data']['filename']

Example JSON:

```
{"filename": ""}
```

IntegralParameters

type object

Required False

JSON Path

- # ['data']['IntegralParameters']

Example JSON:

```
{"IntegralParameters": {}}
```

graphs

Type array() items: {*kind*, *conf*, *columnLabels*, *array*}

Required False

JSON Path

- # ['data']['graphs']

Example JSON:

```
{"graphs": []}
```

kind

Type string

values [Radial, Slice]

Required False

JSON Path

- # ['data']['graphs'][0]['kind']

Example JSON:

```
{"kind": "Radial"}
```

conf

Type object

Required False

JSON Path

- # ['data']['graphs'][0]['conf']

Example JSON:

```
{"conf": null}
```

columnLabels

Type array() items: string

Required False

JSON Path

- # ['data']['graphs'][0]['columnLabels']

Example JSON:

```
{"columnLabels": []}
```

array

Type array() items:

Required False

JSON Path

- # ['data']['graphs'][0]['array']

Example JSON:

```
{"array": []}
```

3.6.3 Leash Datamerge and saxsdmerge Schema

The data merge tool can merge (join) log files with data from the images and create consolidated data files for archiving an further processing.

The '*' signifies a required Field.

Type object

Contains *TimeOffset**, *LogDataTables*, *OutputFormats**, *OutputFileName**, *HDFOptions**

Required False

JSON Path

- #

Example JSON:

```
{
  "TimeOffset": "0",
  "HDFOptions": {
    "IncludeCHI": false,
    "IncludeTIF": false
  },
  "OutputFormats": {
    "exel": false,
    "hdf": false,
    "json": false,
    "csv": false
  },
  "OutputFileName": "../results/merged"
}
```

TimeOffset

If offset is not found otherwise, use this as offset for all log data

Type number in Seconds

Required True

Default 0

JSON Path

- # ['TimeOffset']

Example JSON:

```
{"TimeOffset": "0"}
```

LogDataTables

Define log files to consolidate with image data. If more than one defined, they will be joined and missing values will be interpolated.

Type array() items: {*TimeEpoch*, *TimeOffset*, *FirstImageCorrelation*, *Name*, *Files*}

Required False

JSON Path

- # ['LogDataTables']

Example JSON:

```
{"LogDataTables": []}
```

TimeEpoch

Time epoch

Type string

values [Mac, Unix]

Required True

Default Mac

JSON Path

- # ['LogDataTables'][0]['TimeEpoch']

Example JSON:

```
{"TimeEpoch": "Mac"}
```

TimeOffset

If offset is not found otherwise, use this as offset.

Type number in Seconds

Required True

Default 0

JSON Path

- # ['LogDataTables'][0]['TimeOffset']

Example JSON:

```
{"TimeOffset": "0"}
```

FirstImageCorrelation

Find offset for all log data by correlating first image with first entry of this table.

Type boolean

Required True

Default False

JSON Path

- # ['LogDataTables'][0]['FirstImageCorrelation']

Example JSON:

```
{"FirstImageCorrelation": false}
```

Name

Name field to be used as prefix in the joined column names.

Type string

Required True

Default log

JSON Path

- # ['LogDataTables'][0]['Name']

Example JSON:

```
{"Name": "log"}
```

Files

One log table may be one file, or a list of files to be concatenated.

Type array() items: {*RemotePath*, *LocalPath*}

Required False

JSON Path

- # ['LogDataTables'][0]['Files']

Example JSON:

```
{"Files": []}
```

RemotePath

Path of logfile on server if used in server mode.

Type array() items: string

Required False

JSON Path

- # ['LogDataTables'][0]['Files'][0]['RemotePath']

Example JSON:

```
{"RemotePath": []}
```

LocalPath

Path of logfile on client. Overrides *RemotePath*.

Type string

Required False

JSON Path

- # ['LogDataTables'][0]['Files'][0]['LocalPath']

Example JSON:

```
{"LocalPath": ""}
```

OutputFormats

List of outputformats to write the consolidated log or the consolidated 'hdf' file.

Type object

Contains *csv**, *hdf**, *exel**, *json**

Required True

JSON Path

- # ['OutputFormats']

Example JSON:


```
{
  "OutputFormats": {
    "exel": false,
    "hdf": false,
    "json": false,
    "csv": false
  }
}
```

csv

Type boolean

Required True

Default False

JSON Path

- # ['OutputFormats']['csv']

Example JSON:

```
{"csv": false}
```

hdf

Type boolean

Required True

Default False

JSON Path

- # ['OutputFormats']['hdf']

Example JSON:

```
{"hdf": false}
```

exel

Type boolean

Required True

Default False

JSON Path

- # ['OutputFormats']['exel']

Example JSON:

```
{"exel": false}
```

json

Type boolean

Required True

Default False

JSON Path

- # [*OutputFormats*][*json*]

Example JSON:

```
{"json": false}
```

OutputFileName

Type string

Required True

Default ../results/merged

JSON Path

- # [*OutputFileName*]

Example JSON:

```
{"OutputFileName": "../results/merged"}
```

HDFOptions

Options only relevant to hdf export.

Type object

Contains *IncludeCHI**, *IncludeTIF**

Required True

JSON Path

- # [*HDFOptions*]

Example JSON:

```
{"HDFOptions": {"IncludeCHI": false, "IncludeTIF": false}}
```

IncludeCHI

Whether to include the .chi files as strings.

Type boolean

Required True

Default False

JSON Path

- # [*HDFOptions*][*IncludeCHI*]

Example JSON:

```
{"IncludeCHI": false}
```

IncludeTIF

Whether to include the images as integer array.

Type boolean

Required True

Default False

JSON Path

- # [*HDFOptions*]['*IncludeTIF*']

Example JSON:

```
{"IncludeTIF": false}
```


THE DEECTOR CALIBRATION FILE

The `SAXS.calibration` class and the `saxsdog` tool accepts a input file with the calibration data. This input file is written as a JSON file. This is a common syntax to express structured data as text. You might want to read a bit about it before moving on.

4.1 Calibration file Reference

The '*' signifies a required Field.

The SAXS configuration file specifies the parameters of a SAXS sensor calibration. It is written in the JSON format which governs the general syntax.

Type object

Contains *Title*, *Geometry**, *Masks**, *Slices*, *Wavelength**, *PolarizationCorrection*, *Directory**, *Threads*

Required True

JSON Path

- #

Example JSON:

```
{
  "Geometry": {
    "Tilt": {
      "TiltRotDeg": 0,
      "TiltAngleDeg": 0
    },
    "Imagesize": [
      1000,
      900
    ],
    "BeamCenter": [
      800.0,
      400.0
    ],
    "PixelSizeMicroM": [
      100.0
    ],
    "DedectorDistanceMM": 1000.0
  },
  "Wavelength": 1.54,
  "Directory": [],
  "Masks": []
}
```

4.1.1 Title

Type string

Required False

JSON Path

- # ['Title']

Example JSON:

```
{"Title": ""}
```

4.1.2 Geometry

Type object

Contains *Tilt**, *BeamCenter**, *DedectorDistanceMM**, *PixelSizeMicroM**, *Imagesize**

Required True

JSON Path

- # ['Geometry']

Example JSON:

```
{
  "Geometry": {
    "Tilt": {
      "TiltRotDeg": 0,
      "TiltAngleDeg": 0
    },
    "Imagesize": [
      1000,
      900
    ],
    "BeamCenter": [
      800.0,
      400.0
    ],
    "PixelSizeMicroM": [
      100.0
    ],
    "DedectorDistanceMM": 1000.0
  }
}
```

4.1.3 Tilt

The sensor, usually is not perfectly perpendicular to the ray direction. The tilt angle can be specified by giving the following parameters.

Type object

Contains *TiltRotDeg**, *TiltAngleDeg**

Required True

JSON Path

- # ['Geometry']['Tilt']

Example JSON:

```
{"Tilt": {"TiltRotDeg": 0, "TiltAngleDeg": 0}}
```

4.1.4 TiltRotDeg

This gives the angel of the tilt direction.

Type number in degree

Required True

Default 0

JSON Path

- # ['Geometry']['Tilt']['TiltRotDeg']

Example JSON:

```
{"TiltRotDeg": 0}
```

4.1.5 TiltAngleDeg

This gives the angle between the ray direction and the normal to the sensor plane.

Type number in degree

Required True

Default 0

JSON Path

- # ['Geometry']['Tilt']['TiltAngleDeg']

Example JSON:

```
{"TiltAngleDeg": 0}
```

4.1.6 BeamCenter

Gives the beam center in pixel coorinates.

Type array(2) items: number

Required True

Default [800.0, 400.0]

JSON Path

- # ['Geometry']['BeamCenter']

Example JSON:

```
{"BeamCenter": [800.0, 400.0]}
```

4.1.7 DedectorDistanceMM

Distance between diffraction center and sensor.

Type number in Millimeters

Required True

Default 1000.0

JSON Path

- # ['Geometry']['DedectorDistanceMM']

Example JSON:

```
{"DedectorDistanceMM": 1000.0}
```

4.1.8 PixelSizeMicroM

The pixel size on the sensor.

Type array(2) items: number

Required True

Default [100.0]

JSON Path

- # ['Geometry']['PixelSizeMicroM']

Example JSON:

```
{"PixelSizeMicroM": [100.0]}
```

4.1.9 Imagesize

Size of sensor image in pixel.

Type array(2) items: integer

Required True

Default [1000, 900]

JSON Path

- # ['Geometry']['Imagesize']

Example JSON:

```
{"Imagesize": [1000, 900]}
```

4.1.10 Masks

Type array() items: {*MaskFile*, *Oversampling*, *PixelPerRadialElement*, *Name*, *qStart*, *qStop*}

Required True

JSON Path

- # ['Masks']

Example JSON:

```
{"Masks": []}
```

4.1.11 MaskFile

Path of Maskfile

Type string

Required True

Default AAA_integ.msk

JSON Path

- # ['Masks'][0]['MaskFile']

Example JSON:

```
{"MaskFile": "AAA_integ.msk"}
```

4.1.12 Oversampling

Oversampling factor for radial integration. The higher, the longer the setup but the higher the accuracy. More than 3 is probably overkill.

Type integer

Required True

Default 3

JSON Path

- # ['Masks'][0]['Oversampling']

Example JSON:

```
{"Oversampling": 3}
```

4.1.13 PixelPerRadialElement

Expresses the width of a radial step in terms of pixels. '1' means $\delta R \approx 1 \text{ PixelSizeMicroM}$.

Type number in Pixel

Required True

Default 1

JSON Path

- # ['Masks'][0]['PixelPerRadialElement']

Example JSON:

```
{"PixelPerRadialElement": 1}
```

4.1.14 Name

Name for mask configuration.

Type string

Required False

JSON Path

- # ['Masks'][0]['Name']

Example JSON:

```
{"Name": ""}
```

4.1.15 qStart

Starting q-value for integral parameters.

Type number in nm^{-1}

Required True

Default 0

JSON Path

- # ['Masks'][0]['qStart']

Example JSON:

```
{"qStart": 0}
```

4.1.16 qStop

Ending q-value for integral parameters.

Type number in nm^{-1}

Required True

Default 0

JSON Path

- # ['Masks'][0]['qStop']

Example JSON:

```
{"qStop": 0}
```

4.1.17 Slices

Slices are designed to analyse GISAXS data. It allows you to specify slices along the sensor axis and get intensity along q_x , q_z directions.

Type array() items: {*Direction*, *Plane*, *CutPosition*, *CutMargin*, *IncidentAngle*, *MaskRef*}

Required False

JSON Path

- # ['Slices']

Example JSON:

```
{"Slices": []}
```

4.1.18 Direction

'x' or 'y' direction in sensor pixel coordinates.

Type string

values [x, y]

Required True

Default x

JSON Path

- # ['Slices'][0]['Direction']

Example JSON:

```
{"Direction": "x"}
```

4.1.19 Plane

Whether the direction is in plane with scattering surface or vertical to it.

Type string

values [InPlane, Vertical]

Required True

Default InPlane

JSON Path

- # ['Slices'][0]['Plane']

Example JSON:

```
{"Plane": "InPlane"}
```

4.1.20 CutPosition

Cut position in pixel coordinates in the other coordinate then specified in 'Direction'. Origin is top left Corner.

Type integer in Pixel

Required True

Default 0

JSON Path

- # ['Slices'][0]['CutPosition']

Example JSON:

```
{"CutPosition": "0"}
```

4.1.21 CutMargin

Number of pixels left and right from cut to include into the average .

Type integer in Pixel

Required True

Default 1

JSON Path

- # ['Slices'][0]['CutMargin']

Example JSON:

```
{"CutMargin": 1}
```

4.1.22 IncidentAngle

Angle of incidence in GISAXS setup.

Type number in degree

Required True

Default 0

JSON Path

- # ['Slices'][0]['IncidentAngle']

Example JSON:

```
{"IncidentAngle": 0}
```

4.1.23 MaskRef

Chose which mask to use for the sclice. '-1' means don't use mask

Type integer

Required True

Default 0

JSON Path

- # ['Slices'][0]['MaskRef']

Example JSON:

```
{"MaskRef": 0}
```

4.1.24 Wavelength

Refined wavelength.

Type number in Angstrom

Required True

Default 1.54

JSON Path

- # ['Wavelength']

Example JSON:

```
{"Wavelength": 1.54}
```

4.1.25 PolarizationCorrection

The scattering direction id dependend on the light polarization. This may be accounted for with the polarization correction.

Type object

Contains *Fraction**, *Angle**

Required False

Default OrderedDict([(u'Fraction', 0.95), (u'Angle', 0)])

JSON Path

- # ['PolarizationCorrection']

Example JSON:

```
{"PolarizationCorrection": {"Angle": 0.0, "Fraction": 0.95}}
```

4.1.26 Fraction

Fraction of light polarized in the given (*Angle*) direction.

Type number

Required True

Default 0.95

JSON Path

- # ['PolarizationCorrection']['Fraction']

Example JSON:

```
{"Fraction": 0.95}
```

4.1.27 Angle

Angle of the polarization plane.

Type number in degree

Required True

Default 0.0

JSON Path

- # ['PolarizationCorrection']['Angle']

Example JSON:

```
{"Angle": 0.0}
```

4.1.28 Directory

Directory to take into account for processing images. Given as a list of subdirectories.

Type array() items: string

Required True

JSON Path

- # ['Directory']

Example JSON:

```
{"Directory": []}
```

4.1.29 Threads

Type integer

Required False

Default 2

JSON Path

- # [*Threads*']

Example JSON:

```
{"Threads": 2}
```

THE TECHNOLOGY

5.1 Integration as Matrix-Vector Multiplication

Every SAXS image \mathbf{p} is a list of pixels that have an intensity value. This 2d array might as well be addressed as a vector with all the pixels addressable with one index \mathbf{p}_i .

The integration over pixels that are within a certain radial interval is in any case a weighted sum of some of the pixels.

This weighted sum is a scalar product with another vector containing the weight factors. As only the pixels in a radius interval are counted, most of these factors are 0.

$$r = \mathbf{c} \cdot \mathbf{p}$$

As we intend to do all the radial intervals at once, we write it as a matrix vector product.

$$\mathbf{r} = \mathbf{X} \cdot \mathbf{p}$$

The columns are the weight factors for the i^{th} radial element. Rearranged in the order of the image, this looks like the ring element relevant for the radial Point.

This matrix would be quite big as it has the dimensions $\text{len}(\mathbf{r}) \cdot \text{len}(\mathbf{p})$. Fortunately most of the entries are 0 and we can use a sparse matrix representation which uses only about $\sim \text{len}(\mathbf{p})$ of memory, as every pixel is counted only once, or, as we will see, about once.




Figure 5.1: The vector \mathbf{c} displayed as image.

Figure *CircleNoAA* Shows the data of such a matrix column.

5.2 Over sampling




Figure 5.2: Ring with anti aliasing / over sampling.

A pixel might lie on the border of two radial intervals, making it unclear to which one it should be added. By only choosing the nearest one, one may get artifacts in the resulting curve especially when only few pixels contribute. So, how could we calculate to which fraction a pixel should account to one radial interval?

The idea here is to use an algorithm comparable to anti aliasing in computer graphics. We will divide a much larger picture into the radial intervals and down sample it to the real pixels. Which results in nicely balanced factors for the border pixels that add up nicely over joining intervals such that the intensity is conserved. If one looks closer at image *Ring with anti aliasing / over sampling.*, one sees that the ring has soft edges. Quite as it would have through anti aliasing.

5.3 The Geometry

The plane of the sensor is not perfectly normal to the beam. So in order to calculate which pixel is on which cone in the diffracted light, we need to express the geometry somehow.

Every pixel has the polar coordinates r, ϕ with the projected diffraction center in the origin. For each pixel (P) the triangle S,C,P (Sample, Center, Pixel, θ, β, γ) can be fully expressed with the law of cosines.

l is the distance the light travels from the diffraction center to the sensor.

$$l^2 = d^2 + r^2 - 2dr \cos(\pi/2 + \alpha)$$

r is the radial coordinate of the pixel P.

$$r^2 = l^2 + d^2 - 2ld \cos(\theta)$$

from these two formulas the diffraction angle (here) θ can be computed.

$$\theta = \arccos(-r^2 - l^2 - d^2/2ld)$$

α comes from the following relation in figure *Angle between two planes.*




Figure 5.3: The SCP triangle.




Figure 5.4: Angle between two planes.

The angle between the sensor plane and the normal plane to the ray is given by τ . The slope s derived from τ is

$$s = \sin(\tau)$$

On the (red) unit circle in the plane of the sensor the distance to the plane normal to the ray is expressed as

$$h = \sin(\phi)s$$

The angle α is therefore:

$$\alpha = \arcsin(\sin(\tau)\sin(\phi))$$


in Python code this is `SAXS.calc_theta()`:


```
def calc_theta(r, theta, d, tilt, tiltdir):
    alpha=np.arcsin(np.sin(tilt)*np.sin(theta+tiltdir))
    lsquared=d**2 +r**2 -2*d*r*np.cos(np.pi/2+alpha)
    return np.arccos(-(r**2-lsquared-d**2)/(2*np.sqrt(lsquared)*d))
```

This angle θ then is calculated for every sub pixel in the sensor. This number then can be rescaled and rounded to the nearest integer in order to get unique integer labels for all the pixels. This labels are the index of the radial interval.

5.3.1 Tilt Angle Correction Test

To check if the tilt angle correction is working, lets create some fake calibration data, with the following peaks in the diffraction curve:





This was done with this configuration file:

```
{
  "Geometry": {
    "Tilt": {
      "TiltAngleDeg": -10,
      "TiltRotDeg": 73.569
    },
    "DetectorDistanceMM": 1031.657,
    "BeamCenter": [
      808.37,
      387.772
    ],
    "Imagesize": [
      1043,
      981
    ],
    "PixelSizeMicroM": [
      172.0,
      172.0
    ]
  },
  "Directory": [
    "."
  ],
  "Masks": [
    {
      "PixelPerRadialElement": 1,
      "MaskFile": "emptymask.tif",
      "Oversampling": 2,
      "qStart": 0.0,
      "qStop": 5.0
    }
  ]
}
```

```
    },  
    ],  
    "Wavelength": 1.54  
  }  
}
```

Which amounts to a large tilt, lets see what Fit2d makes of it

```
INFO: SOLUTION 2  
INFO: Best fit beam centre (X/Y mm) = 66.78356 138.9544  
INFO: Best fit beam centre (X/Y pixels) = 388.2765 807.8745  
INFO: Cone 1 best fit 2 theta angle (degrees) = 1.392646  
INFO: Cone 2 best fit 2 theta angle (degrees) = 2.780810  
INFO: Cone 3 best fit 2 theta angle (degrees) = 4.168858  
INFO: Cone 4 best fit 2 theta angle (degrees) = 5.556975  
INFO: Cone 5 best fit 2 theta angle (degrees) = 6.944765  
INFO: Best fit angle of tilt plane rotation (degrees) = 73.59509  
INFO: Best fit angle of tilt (degrees) = -10.00601  
INFO: Estimated coordinate radial position error (mm) = 0.7136476E-02  
INFO: Estimated coordinate radial position error (X pixels) = 0.4149114E-01
```

Seems OK.

5.4 Polarization Correction

The polarization correction is expected to be small at small angles, but it is deemed important.

$$I_{cor} = I_j [P(1 - (\sin(\phi)\sin(2\theta))^2)(1 - P)(1 - (\cos(\phi)\sin(2\theta))^2)]$$


where ϕ is the azimuthal angle on the detector surface (defined here clockwise, 0 at 12 o'clock) 2θ the scattering angle, and P the fraction of incident radiation polarized in the horizontal plane (azimuthal angle of 90°) The polarization correction is configured by two parameters in *PolarizationCorrection*. Its factors are included in the integration matrix (operator).

This input:


```
{  
  "Geometry": {  
    "Tilt": {  
      "TiltAngleDeg": -0.56,  
      "TiltRotDeg": 73.569  
    },  
    "DetectorDistanceMM": 1031.657,  
    "BeamCenter": [  
      808.37,  
      387.772  
    ],  
    "ImageSize": [  
      1043,  
      981  
    ],  
    "PixelSizeMicroM": [  
      172.0,  
      172.0  
    ]  
  },  
  "Directory": [  
    "."  
  ],  
  "PolarizationCorrection": {  
    "Angle": 0,  
    "Fraction": 1  
  }  
}
```

```
},  
"Masks": [  
  {  
    "PixelPerRadialElement": 1,  
    "MaskFile": "./data/AAA_integ.msk",  
    "Oversampling": 2,  
    "qStart": 0.0,  
    "qStop": 5.0  
  }  
],  
"Wavelength": 1.54  
}
```


Gives:



If the correction factors are all correctly in the algorithm, the integration of an image containing $1/I_{corr}$ should give constant 1.0.




Just for checking: integrating a picture with only ones gives something different:




This are the wiggles that come from the polarization correction pattern

5.5 Compare With Fit2d



The program fit2d, which this package aims to partly replace, is the standard, so we better include a comparison plot here:




Okay, this doesn't show much but if we plot the difference:



Still looks okay.



In the tail region the blue halo (Poisson error) signifies that there are not enough counts to make good statistics.



5.6 Integrating a Constant Image With Masked Values


This test shows that nothing wrong happens at mask borders. For those we want to integrate an image that is one everywhere except for the masked regions

We use the following calibration without Polarization correction and mask:


```
{
  "Geometry": {
    "Tilt": {
      "TiltAngleDeg": -0.56,
      "TiltRotDeg": 73.569
    },
    "DetectorDistanceMM": 1031.657,
    "BeamCenter": [
      808.37,
      387.772
    ],
    "Imagesize": [
      1043,
      981
    ],
    "PixelSizeMicroM": [
      172.0,
      172.0
    ]
  },
  "Directory": [
    "."
  ],
}
```

```
"Masks": [  
  {  
    "PixelPerRadialElement": 1,  
    "MaskFile": "data/AAA_integ.msk",  
    "Oversampling": 2,  
    "qStart": 0.0,  
    "qStop": 5.0  
  }  
],  
"Wavelength": 1.54  
}
```

The image we are going to integrate is exactly the array the `SAXS.openmask()` returns:



The result is constant 1 (where the intensity is not 0), save $2e-12$.



Doing the same with Fit2d,


```
$ fit2d -svar\#IN=mask.tif -dim2000x2000 -svar\#OUT=data/const.chi -mac../data/AAA_integ_Pilatus1
```

```
STATUS: The error was identified in module:  
IO: Input/Output and Status
```

```
STATUS: Position where error condition was identified  
Subroutine IO_OPEN_INMACRO V0.2
```

```
STATUS: The error condition has been classified as:  
Failed to open file: (file not found or bad file name ?)  
...
```

results in something similar, just with less precision, about $10e-7$. Probably because of single precision arithmetics.



5.7 Statistics

5.7.1 Poisson Statistics

The most important error is the statistical fluctuation that stems from the randomness of the scattering events. Counts of such events follow the Poisson distribution. Such, the error (σ) is \sqrt{n} for a count of n . The result of which is, that the relative error $\frac{\sqrt{n}}{n}$ rapidly gets small for larger counts.

Each Pixel in the SAXS sensor counts the number of events, and thus follows the Poisson statistics. The error of a sum of pixels is calculated as.

$$\sigma_{sum} = \sqrt{\sum_i \sigma_i^2}$$

which means here

$$\sigma_{sum} = \sqrt{\sum_i n_i}$$


Rescaled over the number of pixels (P) in the sum this gives:

$$\sigma_{sum} = \frac{\sqrt{\sum_{i=1}^P n_i}}{P}$$


The `SAXS.calibration.plot()` method of the `SAXS.calibration` class will give you the Poisson error along with the standard deviation. So for regions, where the total number of counts is too small, you can see if there is a significant error. This might occur, if too few pixels are used for a data point or the intensity is just too small.

5.7.2 Standard Deviation

The standard deviation of the mean that is taken through the integration is not as such particularly useful to estimate the error of the resulting intensities because there are quite a few things that produce an angle dependence. In an optimal case, if the angle dependence can be corrected with the Polarization correction, the standard deviation of the integration might be very small. In an ordinary case the standard deviation gives you a measure of how spread the intensities within a radius interval are.



The standard deviation is bright yellow and the Poisson error is blueisch



If the calibration is wrong you will for example see in the standard deviation. Like in this example. Here the beam center is wrong.

5.8 GISAXS Slices

In grating incident SAXS looking at radial integration doesn't make much sense. In GISAXS you rather want to look at horizontal or vertical slices. That is horizontal or vertical with respect to the scattering surface.

The SAXSDog tools allow to specify *Slices* of pixels and allow to look at them at the q_y , q_z scale as they are used in GISAX analysis.

The `slice` class implements this functionality very much the same as the radial integration except that the labels of the pixels simply are the integer x coordinate in pixels when we want a slice in x direction and the y coordinate when we want the slice in y direction. Oversampling doesn't make any sense in this scenario but the rest is the same. The only subtle thing is to calculate the q_y and q_z scale. Because the detector coordinate system may have the x or the y axis aligned with the scattering surface the *Slices* allows to choose whether we are in plane or vertical to the scattering surface. If you want to get correct q_z values you must also specify the incident angle α_i , (*IncidentAngle*).

SAXSDOG DEVELOPER DOCUMENTATION

6.1 The Leash

The Leash GUI is part of the SAXSdog Network and the most user facing software. Its design allows to extend it with functionality by touching only small well defined parts. One is the Schemas of the configuration files and the other is the `QItemDelegate`. The following sections give an overview of the important parts.

6.1.1 Using the JSON Schema to Extend Leash

This software relies a lot on structured configuration files that have to constantly be checked for validity. This is done by defining the grammar in JSON Schema. This is a language, in its self expressed in JSON, made to specify what values may occur where in the file. This does not only allow for automatically generating documentation, as it is used in this document many times, but you can also use it to generate an GUI that can edit this structured data files.

The tree view in the “Calib” tab of Leash is build by recursively, going to the schema and the data, building the model that can be displayed in the `QTreeView` widget.

So, in order to add new parameters to the view, the only thing you must do, is to add the description to the schema. If you use similar constructs as in the rest of the data, it will work just so.

The leash uses the scheme in `SAXS/schema.json` to build the “Calib” tab and `SAXS/DataConsolidationConf.json` to build the Consolidate tab.

Consider this excerpt of the `SAXS/schema.json`:

```
{
  "type": "object",
  "$schema": "http://json-schema.org/draft-03/schema",
  "required": true,
  "description": "The SAXS configuration file specifies the parameters of a SAXS sensor calibration",
  "additionalProperties": false,
  "properties":
  {
    "Directory":
    {
      "description": "Directory to take into account for processing images. Given as a list of s",
      "type": "array",
      "required": true,
      "minItems": 1,
      "items":
      {
        "type": "string",
        "default": ".",
        "appinfo":
        {
          "editor": "RemoteDir"
        }
      }
    }
  }
}
```

```
    }  
  }  
}
```

The type in the root says `object` which means it is a dictionary like data structure. The possible keys are declared in the `properties` dictionary. Inside the `properties` other nested types are declared in this case an array called `Directory`.

The `appinfo` section is not part of the schema specification it is rather a custom field to tell the application about possible special treatments. In this case we want the Leash use editor widgets to pick a remote directory. This editor widget will be provided by the `QItemDelegate`. Which shall be explained in the following section.

6.1.2 The QItemDelegate

The `QtGui.QTreeView()` class allows to set an item delegate. This happens for example in `calibeditor.py`

```
self.treeview.setItemDelegateForColumn(1,calibeditdelegate.calibEditDelegate( app ))
```

The `calibeditdelegate.calibEditDelegate` class in turn is a custom class derived from `QtGui.QItemDelegate`. This is implemented in `calibeditdelegate.py`.

The constructor is initializing the base class:

```
class calibEditDelegate(QtGui.QItemDelegate):  
    def __init__(self,app, parent=None):  
        super(calibEditDelegate, self).__init__(parent)  
        self.app=app
```

The part that is interesting is the reimplementations of the `createEditor`, `commitAndCloseEditor`, `setEditorData` and `setModelData` methods.

The `createEditor` method is called when the user double clicks an item content cell in the tree view. The default behavior is to make the text content editable but you can return any widget you like, depending on the context.

```
def createEditor(self, parent, option, index):  
    """  
    special method of QItemDelegate class  
    """  
    try:  
        subschema=json.loads(unicode(index.model().data(index,role=im.SUBSCHEMA).toString()))  
    except ValueError:  
        return None  
  
    type= unicode(index.model().data(index,role=im.TYPE).toString())  
    editablearray= unicode(index.model().data(index,role=im.ISEDITABLEARRAY).toString())  
    editortype=None  
    if subschema.get("appinfo"):  
        editortype= subschema.get("appinfo").get("editor")  
  
    print type  
    if "enum" in subschema:  
        isenum="true"  
        enum=subschema['enum']  
    else:  
        isenum="false"  
  
    if type == "integer":  
        spinbox = QtGui.QSpinBox(parent)  
        spinbox.setRange(-200000, 200000)  
        spinbox.setSingleStep(1)
```

```

        spinbox.setAlignment (QtCore.Qt.AlignRight | QtCore.Qt.AlignVCenter)
    return spinbox
elif type == "number":
    spinbox = QtGui.QDoubleSpinBox (parent)
    spinbox.setRange (-200000, 200000)
    spinbox.setSingleStep (0.1)
    spinbox.setDecimals (4)
    spinbox.setAlignment (QtCore.Qt.AlignRight | QtCore.Qt.AlignVCenter)
    return spinbox
elif editablearray=="editablearray":
    arrayeditdialog=arrayeditdialog (index, parent)
    return arrayeditdialog
elif type== "object" or type=="array" or type=="arrayitem" :
    return None
elif isenum=="true":
    combobox = QtGui.QComboBox (parent)
    combobox.addItem (sorted ( enum))
    return combobox
elif editortype=="File":
    dirname= os.path.dirname (unicode (index.model ().filename))
    filepicker=QtGui.QFileDialog (directory=dirname)
    filepicker.setMinimumSize (800,500)
    filepicker.setFileMode (filepicker.ExistingFile)
    return filepicker
elif editortype=="RemoteDir":
    return RemoteDirPicker (self.app, parent, index)
elif editortype=="RemoteFile":
    return RemoteDirPicker (self.app, parent, index, showfiles=True)
else:
    return QtGui.QItemDelegate.createEditor (self, parent, option,
                                             index)

```

The `createEditor` is called with an `index` object. Which is a class that is used by the `QtGui.QStandardItemModel` class to represent the data in a form the tree view can display it.

In the implementation in `jsonschematree.py` the items have the subschema describing themselves and their children stored in special data attributes. We can use this to chose which editor to present to the user, depending of the type and role of the item on hand. Integers get a `QtGui.QSpinBox`, Enumerations get `QtGui.QComboBox` to select one of the options.

In case the item has `File` in the `appinfo/editor` field,

```

{
  "MaskFile":
  {
    "description": "Path of Maskfile",
    "type": "string",
    "default": "AAA_integ.msk",
    "required": true,
    "appinfo":
    {
      "editor": "File",
      "display": "MaskFile"
    }
  }
}

```

double clicking the cell will give the user a file system dialog to select a local file. This time we did not get a small widget that fits into the cell, we got a separate dialog. This means it is possible to launch any kind of fancy dialog from here. Think “mask editor”, “powder diffraction calibration” anything you like.

The `"display": "MaskFile"` field will cause another method to execute custom behavior. The `setModelData` method. In this case it will load the mask file and display the picture in another cell in the tree view.

6.2 The Image Queue

The `imagequeue` class manages how and when to integrate images. It is instantiated by the server when you load up a new calibration and start a new queue. Or, alternatively the `saxsdog` command line tool will also create an image queue. It takes as argument a list of integration recipes e.g. radial integration or slices. These recipes can be any Python object that knows how to do something with images as long they implement the `integratechi()` method with the same API as the others.

In the initialization process, it will create a queue object, which is a very powerful synchronized data structure which can even be accessed by subprocesses. If the server creates it it will also create a process to listen to the Feeder service and push the image paths into the queue of the `imagequeue` object.

For the work to begin the `imagequeue.start()` method needs to be called. This will create the worker subprocesses to consume the images from the queue.

```
for threadid in range(1, self.options.threads):
    print "start proc [" , threadid, "]"
    worker=Process(target=funcworker, args=(self, threadid))
    worker.daemon=True
    self.pool.append(worker)
    worker.start()
```

The `imagequeue` will launch and manage as many workers as configured in the calibration. The workers are in an infinite loop where they wait until a new image arrives through the queue to decide whether they are configured to work on the directory the images are in. If so they will process the image and push a small report into the history queue. This report includes the time (for the histogram) and the files written.

If the `readdir` command is issued to the server, it will call the `imagequeue.fillqueuewithexistingfiles()` method which will fill the queue with all ".tif" files it finds in the configured directory.

6.3 Data Consolidation or Datamerger

After the images are integrated and the measurements are done there remain a few data consolidation tasks. The `datamerge` module provides this functionality. The `Saxsdmerge` commandline tool is one way to use it, the other is via the `Leash`. Its main goal is to merge the logfiles with the parameters logged in the images and detector logs.

The result is a table where each image has one row and the columns are all the available parameters from continuous logs or logs that log not regularly but only if an image is requested.

This is where future versions may include a complete HDF export of all relevant data.

SAXS MODULE API

class `SAXS.calibration` (*config, mask=None, attachment=None*)

This class represents a calibration for SAXS diffraction. After initialization, the `integrate()` method can compute the radial intensity very fast.

Parameters

- **config** (*dict*) – the calibration json object as dictionary *Calibration file Reference*:
- **mask** (*dict*) – The mask of the list to use for this instance
- **attachment** (*dict*) – If mask is not to be read from filesystem it must be provided as attachment

integParameters (*Intensity*)

function to calculate the integral parameters between `qStart` and `qStop` returns `I0, I1, I2`

integrate (*image*)

Integrate a picture.

Parameters `image` (*numpy.array(dim=2)*) – Sensor image to integrate as 2d *NumPy* array

Returns Returns Angle and intensity vector as a tuple (angle,intensity)

integratechi (*image, path, picture*)

Integrate and save to file in “chi” format.

Parameters

- **image** (*np.array()*) – Image to integrate as numpy array
- **path** (*string*) – Path to save the file to

Returns Scattering curve data as numpy array

integrateerror (*image*)

Integrates an image and computes error estimates.

Parameters `image` (*np.array()*) – Image to integrate as numpy array

Returns The intensity the standard deviation and the Poisson statistics error in a numpy array.

plot (*image, outputfile='', startplotat=0, fig=None*)

Plot integrated function for image in argument.

Parameters

- **image** (*numpy.array(dim=2)*) – Sensor image to integrate as 2d *NumPy* array
- **outputfile** (*string*) – File to write plot to. Might be any image format supported by matplotlib.
- **startplotat** (*integer*) – radial point from which to start the plot

polcorr (*Pfrac, rot*)

Polarization Correction

class SAXS.**slice** (*conf, sliceconf, attachments=[]*)
Initializes slice integration

Parameters

- **conf** (*dictionary*) – Detector calibration
- **sliceconf** (*dictionary*) – The part of detector calibration which deals with this slice.
- **attachments** (*dictionary*) – Mask data.

integratechi (*image, path, picture*)
Integrate and save to file in “chi” format.

Parameters

- **image** (*np.array()*) – Image to integrate as numpy array
- **path** (*string*) – Path to save the file to

Returns Scattering curve data as json structure

plot (*image, outputfile='', startplotat=0, fig=None*)
dummy function in order to not trip up image queue

SAXS.**calc_theta** (*r, phi, d, tilt, tiltdir*)

Calculates the diffraction angle from pixel coordinates. It does work when called with arrays. See [Geometry](#)

Parameters

- **r** (*float*) – Distance to beamcenter.
- **phi** (*float*) – Angle[rad] from polar sensor plane coordinates.
- **d** (*float*) – distance to diffraction center.
- **tilt** (*float*) – Angle[rad] of sensor plane tilt.
- **tiltdir** (*float*) – Angle[rad] of direction of tilt.

Returns theta

SAXS.**scalemat** (*Xsize, Ysize, ov*)

Computes a scaling projection for use in computing the pixel weights for integration

Parameters

- **Xsize** (*int*) – Picture size in X direction.
- **Ysize** (*int*) – Picture size in Ydirection.
- **ov** (*int*) – Number of oversampling ticks in x ynd y direction
- **corr** (*array*) – Polarizationn an other correction factors

Returns sparce matrix toing the scaling

SAXS.**openmask** (*mfile, attachment=None*)

Open the mask file especially the *.msk file. Unfortunately there is no library module for msk files available also no documentation. So, for the msk file, we have a very brittle hack it works for our sensor. Nevermind any other resolution or size.

Parameters **config** (*object*) – Calibration config object.

Returns Mask as logical numpy array.

SAXS.**convert** ()

This implements the functionality oft [The Converter](#). It parses the commandline options and converts the Fit2d info file to the JSON data used by the SAXS.callibration class.

SAXS.**saxsdog** ()

This implements the functionality of [The Saxsdog](#)

class SAXS.**imagequeue** (*Cals, options, directory, conf*)

This class keeps a queue of images which may be worked on in threads.

Parameters

- **Cal** (*SAXS.calibration*) – The SAXS Calibration to use for the processing
- **options** (*optparser*) – The object with the comandline options of the saxsdog
- **args** (*list*) – List of command line options

fillqueuewithexistingfiles ()

Fill the queue with the list of images that is already there.

start ()

Start threads and directory observer.

class SAXS.**Server** (*conf, serverid, stopflag=None, serverdir=None*)

class to manage a saxsdog server

authenticate (*data*)

check signature of request

commandhandler (*object, attachment*)

start ()

start server loop

start_image_queue (*object, attachment*)

prepare new image queue start processing threads

SAXS.**initcommand** (*options, arg, conf*)

Interface for issuing leash commands

SAXS.**validateResponse** (*message*)

Validate response from saxsdog server against the schema.

Get this guide as PDF.



A

authenticate() (SAXS.Server method), 67

C

calc_theta() (in module SAXS), 66

calibration (class in SAXS), 65

commandhandler() (SAXS.Server method), 67

convert() (in module SAXS), 66

F

fillqueuewithexistingfiles() (SAXS.imagequeue method), 67

I

imagequeue (class in SAXS), 66

initcommand() (in module SAXS), 67

integParameters() (SAXS.calibration method), 65

integrate() (SAXS.calibration method), 65

integratechi() (SAXS.calibration method), 65

integratechi() (SAXS.slice method), 66

integrateerror() (SAXS.calibration method), 65

O

openmask() (in module SAXS), 66

P

plot() (SAXS.calibration method), 65

plot() (SAXS.slice method), 66

polcorr() (SAXS.calibration method), 65

S

SAXS (module), 65

saxsdog() (in module SAXS), 66

scalemat() (in module SAXS), 66

Server (class in SAXS), 67

slice (class in SAXS), 65

start() (SAXS.imagequeue method), 67

start() (SAXS.Server method), 67

start_image_queue() (SAXS.Server method), 67

V

validateResponse() (in module SAXS), 67